




TYC12-TYC14 >>>

CODE	PRODUCT	DESCRIPTION	PACKING
TYC12	FILTER TIP 10ul WB	With Box	96/rack, 10racks/unit, 5units/ctn
TYC13	FILTER TIP 200ul WB	With Box	96/rack, 10racks/unit, 5units/ctn



TYJ07 >>>

CODE	PRODUCT	DESCRIPTION	PACKING
TYJ07	BLOOD LANCET SAFETY PA-SS	Pressure Activated, Stainless Steel, Sterile (gamma radiation) Gauge 23g	100pcs/box, 20box/ctn



TYJ11 >>>

CODE	PRODUCT	DESCRIPTION	PACKING
TYJ11	FEMALE SWAB 3cm	E.O/Gamma Sterilization. Stick plastic.	100pcs/bag, 200box/ctn
TYJ11	FEMALE SWAB 8cm	E.O/Gamma Sterilization. Stick plastic.	100pcs/bag, 200box/ctn



GOLD CAP SST TUBE >>>

GOLD CAP tube contains a clot activator and serum gel separator. This tube can be used for various chemistry, serology and immunology tests. After centrifugal, the specimen will be separated into three parts, namely: serum, separation gel and blood corpuscle. The thick layer of separation gel ensures complete isolation of the serum from the blood corpuscle.

TYK04-TYK06 >>>

CODE	MATERIAL	ADDITIVE	SIZE (mm)	VOL (ml)	PACKING	CTN SIZE (cm)	GW (kg)
TYK05	PET	Clot Activator & Serum Gel Separator	13*100	5	100pcs/rack 1800pcs/ctn	52*38*38	14.4
TYK06	PET	Clot Activator & Serum Gel Separator	16*100	8.5	100pcs/rack 1200pcs/ctn	51*41*25	18



BLUE CAP PT TUBE >>>

BLUE CAP PT tube contains a sodium citrate solution, suitable to routine clinical examination. The concentration of sodium citrate is 109mmol/L and the accurate blood to anticoagulant ration is 1:9

TYK10-TYK10 >>>

CODE	MATERIAL	ADDITIVE	SIZE (mm)	VOL (ml)	PACKING	CTN SIZE (cm)	GW (kg)
TYK10	PET	Sodium Citrate 1:9	13*75	1.8	100pcs/rack 1800pcs/ctn	52*38*38	12.6



PURPLE CAP EDTA TUBE >>>

PURPLE CAP tube contains EDTA (Ethylenediametraacetic Acid) as an anticoagulant. This tube can be utilized for most haematological procedures and are the preferred method for molecular tests.

TYK12-TYK13 >>>

CODE	MATERIAL	ADDITIVE	SIZE (mm)	VOL (ml)	PACKING	CTN SIZE (cm)	GW (kg)
TYK12	PET	EDTA K2/K3	13*100	5	100pcs/rack 1800pcs/ctn	52*38*38	12.4
TYK13	PET	EDTA K2/K3	16*100	8.5	100pcs/rack 1200pcs/ctn	51*41*25	14



GREEN CAP HEPARIN TUBE >>>

GREEN CAP tube contains lithium heparin and a gel separator used for the collection of heparinized plasma for routine and emergency chemistry tests.

TYK15-TYK16 >>>

CODE	MATERIAL	ADDITIVE	SIZE (mm)	VOL (ml)	PACKING	CTN SIZE (cm)	GW (kg)
TYK15	PET	HEPARIN LITHIUM	13*100	5	100pcs/rack 1800pcs/ctn	52*38*38	14.4
TYK16	PET	HEPARIN LITHIUM	16*100	8.5	100pcs/rack 1200pcs/ctn	51*41*25	14



GREY CAP GLUCOSE TUBE >>>

GREY CAP tube contains potassium oxalate and sodium fluoride as a preservative. This tube is used to preserve glucose in whole blood and the original characters of the blood specimen can be maintained for an extended period.

TYK18 >>>

CODE	MATERIAL	ADDITIVE	SIZE (mm)	VOL (ml)	PACKING	CTN SIZE (cm)	GW (kg)
TYK18	PET	Sodium Fluoride & Potassium Oxalate	13*100	5	100pcs/rack 1800pcs/ctn	52*38*38	15.8

